

## anti-human CD97

The antibody MEM-180 recognizes an unique epitope on CD97, a 75-85 kDa surface glycoprotein of G-protein-coupled receptor family, expressed on activated B and T lymphocytes, monocytes/macrophages, dendritic cells and granulocytes.

**CD97** is a G-protein-coupled seven-span transmembrane adhesive receptor that is constitutively expressed on granulocytes and monocytes and rapidly upregulated on T and B cells upon activation. CD97 is produced in alternatively spliced forms and its cellular ligand is CD55 (DAF), which protects various cell types from complement-mediated damage. Interaction of CD97 on leukocytes and CD55 on vessel cells probably facilitate leukocyte activation and migration into the tissues, similarly, CD97 seems to play a role in tumour migration and invasiveness. CD97 is involved in T cell regulation and peripheral granulocyte homeostasis.

<b>Clone:</b>	MEM-180
<b>Isotype:</b>	Mouse IgG1
<b>Physical state:</b>	Purified from ascites by protein-A affinity chromatography
<b>Purity:</b>	> 95% (by SDS-PAGE)
<b>Buffer:</b>	PBS with 15 mM sodium azide (pH 7.4)
<b>Storage Instruction:</b>	Store at 4 °C. For long-term storage aliquot and store at -20°C. Avoid freeze/thaw cycles.
<b>Application:</b>	Flow Cytometry, Immunoprecipitation
<b>Reference:</b>	Leukocyte Typing VI., Kishimoto T. et al. (Eds.), Garland Publishing Inc. (1997)
<b>Quantity:</b>	0.1 mg
<b>Order N°:</b>	H12211

**Warning:** Sodium azide is harmful if swallowed (R22). Keep out of reach of children (S2). Keep away from food, drink and animal feeding stuff (S13). Wear suitable protective clothing (S36). If swallowed, seek medical advice immediately and show this container or label (S46). Contact with acids liberates very toxic gas (R32).

This material is offered for **research only**. Not for use in human. For in vitro use only. EuroBioSciences will not be held responsible for patent infringement or other violations that may occur with the use of our products.