



[Home](#)

Recommendations

please note:

The amount of antibodies you need for your assay is user-dependent: (dependent from cell amount, receptor density, incubation volume, incubation time, temperature, quality of the antibodies, quality of the laser and detectors).

You can do between 200 to 500 experiments with each antibody.

If you use the antibodies at the first time, please do some dilutions to find out the best concentration for your individual requirements. We recommend to start with 2µl to 5µl per 1 million cells.

We advise to use our antibodies 4-5µl of the conjugated antibody to 1 million cells

[back to the website](#)